First results of the JUPITER Registry on long-term performance and safety of the Transapical JenaValve

euro

Stephan Ensminger, MD PhD Heart and Diabetes Center NRW, Ruhr-University Bochum, Bad Oeynhausen

On behalf of the JUPITER registry investigators

Utz Kappert, Ardawan Rastan, Ulrich Franke, Hendrik Treede, Walter Eichinger, Florian Rüter, Rüdiger Lange, Friedhelm Beyersdorf, Thom De Kroon, Olaf Wendler









- TA-JenaValve received CE-Mark in September 2011
- Post-market registry to evaluate 5-year long-term safety and effectiveness in 180 elderly patients
- Enrollment of consecutive patients to avoid selection bias "Real world" population according to IFU
- Acute and 30-day interim results for first half of the patients

PCR 2013



Study design

Methodology

- Intent-to-treat analysis
- Adjudication of major VARC events by independent medical reviewer
- 100% SAE monitoring
- Endpoint definition according to VARC

Endpoints

- Primary 30-day mortality
- Secondary Safety

Device success Effectiveness Quality of life (SF-12)

PCR Image: Constraint of the second seco

Access route Deployment Stent Material Valve Material Skirt Material # of Valve Sizes Annulus Range **Delivery System** Features

Transapical Self Expanding Nitinol Native Porcine Aortic Valve Porcine Pericardium 3 (23, 25, 27) 21-27 mm Cathlete Catheter (sheathless), 3 Step Mechanism

- Feeler guided, anatomically correct positioning
- JenaValve Clipping Mechanism
- No rapid pacing needed during prosthesis implantation





3 Step Delivery

True anatomically correct positioning with predefined landing zone



Step 1

- Unsheating of feelers
- Orientate the markers
- Tactile feedback
- Anatomically correct positioning

Step 2

- Clipping the JenaValve on the native leaflets
- Less material in the LVOT
- Free left and right ostium

Step 3

- Final valve release
- Correct positioning
- No rapid pacing needed during procedure





Baseline Characteristics	n=88
Age (years)	80.8 ± 6.1
Female	38.6 % (34)
NYHA class ≥III	85.2 % (75)
Coronary artery disease	46.6 % (41)
Previous CABG	15.9 % (14)
Chronic renal insufficiency	39.8 % (35)
LVEF (%)	51.5 ± 11.6
Logistic EuroSCORE (%)	24.9 ± 13.5
STS Score	6.1 ± 4.1





Procedural Outcomes

Procedural Outcomes	n=88
Native annulus size (mean \pm SD)	23.6 ± 1.9
Valve sizes: 23	14.8 % (13)
25	47.7 % (42)
27	37.5 % (33)
Duration of TAVI procedure (min)*	8.5 ± 6.6
Procedural success	95.5 %(84)
Valve in valve	2.3 % (2)
Conversion to AVR	2.3 % (2)
Pacemaker implantation (new onset)	12.5 % (11)

* Time from insertion until removal of delivery system









30 day Clinical Outcomes (according to VARC)

Major events	% (n)
Major stroke	0.0% (0)
Spontaneous myocardial infarction	1.3 % (1)
New onset acute kidney injury (<72h, stage 3)	2.7% (2)
Repeat procedure for valve related dysfunction*	1.3 % (1)

* AVR due to moderate PVL





Echocardiography



PCR 2013

Paravalvular Leakage

Post Market Registry to evaluate long-term safety and effectiveness of JenaValve (2nd generation TAVI device)

Unique delivery/fixation system of the JenaValve results in:

- High procedural success
- No major stroke
- Excellent haemodynamics
- Low incidence of paravalvular leakage 97.6 % ≤ mild; 0 % severe